

VET CHROMA

POINT OF CARE PARA IMUNOENSAIOS.



HORMÔNIOS	T4	Tiroxina (T4) Canina e Felina
	CORTISOL	Cortisol Canino e Felino
	cTSH	Hormônio Estimulante de Tireoide (TSH) Canino
	PROGESTERONA	Progesterona Canino
INFLAMATÓRIOS	c-CPR	Proteína C-Reativa Canina
	fSAA	Soro Amiloide A Felino
	cPL	Lipase Pancreática Canina
	fPL	Lipase Pancreática Felina
CARDÍACOS	D-Dímero	D-Dímero canino
	cNT-ProBNP	NT-ProBNP Canino
	fNT-ProBNP	NT-ProBNP Felino
	cTN-I	Troponina I Canina
HEPÁTICO	AFP	Alfa-fetoproteína
RENAL	SDMA	Dimetilarginina simétrica


Deteção e quantificação de biomarcadores essenciais no diagnóstico e acompanhamento de doenças **cardíacas, inflamatórias, hormonais e muito mais!**

Procedimento
geral dos testes
no Vet Chroma:




3 minutos

Resultados Quantitativos
a partir de 3 min.



Técnica fluorescente de alta sensibilidade

- ✓ Alta precisão e reprodutibilidade
- ✓ Resultados entre 3-15 minutos
- ✓ Método de teste fácil e simples
- ✓ Nenhum custo de manutenção
- ✓ Amplo menu de testes em um único equipamento portátil



Step 1



Step 2

16-Nov-2020 / 09:27:07

Step 2

Owner

Test: Progesterone

Lot / Exp: PGQHB03V / 2022.04.11

Comments

Animal


Animal ID

Species: Dog

Dog

Home Next

Step 3



Result

16-Nov-2020 / 09:27:47

Progesterone Result

Owner	
Animal	
Animal ID	
Species	Dog
Lot / Exp	PGQHB03V / 2022.04.11
Comments	
Result	> 25.00 ng/mL

Save Save & Print



Acesse aqui
nosso Guia de
Coleta de Sangue

